

SUPERUSER GUIDE: When nurses are ready to be checked off on Baxter competency, go through each step below, asking the nurse to complete each function.

1. "TURN THE PUMP ON" The nurse will:

- Ensure pump is plugged in before turning on. **SUPERUSER NOTE:** If not plugged in, battery will be damaged
- Turn on with ON/OFF key / Perform speaker test
- There are three things to check before pump goes into main screen in 10 seconds
 - Verify battery charge level: battery charge icon displayed in upper left part of screen
 - Verify PERSONALITY status is correct
 - If new patient: Clear patient history by pressing NEW PATIENT and DONE
- To change the personality if the pump is already in main screen: it must be turned off then on again

2. "LOAD THE ADMINISTRATION SET INTO THE PUMP" The nurse will:

- Once administration set is prepared: Close dark blue key clamp
- If using triple channel pump: Press CHANNEL SELECT.
- Press OPEN key. Load Administration Set. Open regulating roller clamp. Ensure no drops are falling into drip chamber

3. "PROGRAM A PRIMARY INFUSION IN RATE-VOLUME: i.e. RATE: 1 ml per hour. VOLUME 50 mL" The nurse will:

- Single channel pump: Press PRIMARY. Triple channel pump: Press CHANNEL SELECT
- Enter RATE then scroll down to VOLUME to be infused. Enter volume
- Press CONFIRM PRIMARY then START

4. "CHANGE THE RATE TO 2 mL per hour" The nurse will:

- Single channel pump: Press PRIMARY Triple channel pump: Press CHANNEL SELECT.
- To change the rate: Do not need to STOP infusion. Scroll to highlight Rate. Enter new rate. Press CONFIRM
- Press START

5. "CHANGE THE CONTAINER: NEW VOLUME 100 mL" The nurse will:

- Single channel pump: Press PRIMARY. Triple channel pump: Press CHANNEL SELECT
- Press STOP.
- Change container.
- Scroll to Volume. Enter new volume.
- Press CONFIRM
- Press START.

SUPERUSER NOTE: For competency checklist, nurse does not have to actually the change container – just change volume on pump

6. "PROGRAM A SECONDARY INFUSION IN RATE-VOLUME: i.e. An antibiotic in 250 mL over 1 hour"

Ask: "What tubing would you use for a secondary infusion?" (Answer: secondary tubing) The nurse will:

- Single channel pump: Press SECONDARY. Triple channel pump: Press CHANNEL SELECT then SECONDARY
- Enter rate and volume. (should be Rate: 250 ml Volume: 250 ml)
- Enable secondary Callback Feature
- Press CONFIRM SECONDARY
- Press START
- Open ON / OFF clamp on secondary set. Verify drops are falling into secondary set drip chamber

SUPERUSER NOTE: When going through the pump functions with a nurse always clear secondary settings before going on. i.e. when in secondary screen press CLEAR for all areas. If you don't, the secondary program will keep popping up

7. CHANGE MODE - Changes the mode of delivery.

“There are three modes in CHANGE MODE that we need to know: ♦LABEL LINE ♦STANDBY ♦ PRIMARY OR SECONDARY VOLUME-TIME” (critical care may need to know other modes as well)

a. ***“LABEL THE LINE BLOOD or TPN” The nurse will:***

- Single channel pump: Press PRIMARY or SECONDARY. Triple channel pump: Press CHANNEL SELECT
- Press STOP then select CHANGE MODE
- Scroll up to select LABEL LINE
- Scroll down to select desired label with arrow key. Press SELECT. Press START. The chosen label will show on the screen

SUPERUSER NOTE: Nurses need to know that they should not use label line for medication names

b. ***“PLACE PUMP IN STANDBY MODE” The nurse will:***

- Single channel pump: Press PRIMARY or SECONDARY. Triple channel pump: Press CHANNEL SELECT
- Press STOP then select CHANGE MODE
- Scroll up to STANDBY
- Press STANDBY. Message will appear: “Do you want to put this channel on standby” Select ‘YES’

SUPERUSER NOTE: Shortcut to standby on triple lumen pumps: double click on channel select button when infusion is stopped

c. ***“SET PUMP TO RUN BLOOD IN PRIMARY VOLUME-TIME: 326 mL OVER 3 HOURS” The nurse will:***

- Single channel pump: Press PRIMARY or SECONDARY. Triple channel pump: Press CHANNEL SELECT
- Press STOP
- Select CHANGE MODE
- Scroll to highlight PRIMARY or SECONDARY VOLUME-TIME. Press SELECT. Enter volume amount and the time period for the infusion in hours / minutes.
- Press CONFIRM PRIMARY or SECONDARY. Rate will be entered automatically.
- Press START

SUPERUSER NOTE: you can ask nurse to program another VOLUME-TIME: ***“Morphine 1 mg/hour”***

8. GUARDIAN LIBRARY (GL)

“PHYSICIAN ORDER: insulin infusion 1 unit per hour: how would you reconstitute”: (answer: GL monograph)

Ask “What tubing would you use?” (Answer: a primary tubing)

SUPERUSER NOTE: You are looking to see that the nurse chooses to use the GL monograph to reconstitute and that they program the pump in GL

SUPERUSER NOTE: Nurses need to know that intermittent infusions are never programmed in gl – only continuous infusions

The nurse will:

- Single channel pump: Press PRIMARY. Triple channel pump: Press CHANNEL SELECT
- Press STOP
- Press COLLEAGUE GUARDIAN. Scroll to highlight desired drug with arrow key
- Press SELECT. All areas will be entered automatically except DOSE
- Compare dosage/concentration on pump with *Colleague Guardian Monograph* to ensure they are identical.
- Enter Dose. Press CONFIRM SETTINGS. Press START
- Guardian limits may be viewed in this window by pressing VIEW LIMITS
- Symbol for Guardian Feature: Mortar / Pestle in top right corner. Name/ Dose is in yellow lettering.
- An Independent Double Check (2 nurses) required for high alert medications. Both nurses sign the MAR.

SUPERUSER NOTE ask nurse to program another drug in GL if you choose: *i.e.* "Lasix 5 mg per hour"

9. CHANGING THE DOSE IN GUARDIAN

"Change the dose to ___" The nurse will:

- **Single channel pump:** Press PRIMARY. **Triple channel pump:** Press CHANNEL SELECT
- Enter new dose
- Press CONFIRM SETTINGS t
- Press START

10. ENTER A DOSE THAT OVERRIDES GUARDIAN DOSE

"Program Insulin 12 units per hour" The nurse will:

- **Single channel pump:** Press PRIMARY. **Triple channel pump:** Press CHANNEL SELECT
- Enter dose that is outside guardian limits
- Press CONFIRM SETTINGS
- 'LIMITS WARNING' will appear top left corner. The word 'High' or 'Low' will appear beside dose
- Choose 'Accept Dose' or 'Cancel Dose'. If accepting dose: **An Independent Double Check (2 nurses) required to override Guardian limits. Both nurses must sign MAR**
- Press CONFIRM SETTINGS
- Press START. While in override, the dose will be in red lettering

SUPERUSER NOTE : Nurse needs to know that before accepting a dose out of limits, two nurses must verify that the dose is correct

11. ALARMS AND ALERTS:

- Nurse knows how to identify and correct: i.e. ♦ AIR DETECTED ♦ UPSTREAM AND DOWNSTREAM OCCLUSION ALARMS. **THE NURSE CAN REVIEW THIS IN THE OPERATORS MANUAL ON THE DESKTOP**

12. FLOW CHECK GRAPHIC **just show this to the nurse**

- Displays downstream in-line resistance to flow

13. **"UNLOAD ADMINISTRATION SET" The nurse will:**

- Close regulating roller clamp
- Press OPEN
- Remove administration set

14. **"TURN PUMP OFF". The nurse will:**

- Press ON / OFF Key. Press ON / OFF key a second time