

MEMO | NOTE



Staff Development/Infection Prevention & Control

Centre Albert-Galliot

Box 130, 44 Rogers Street

Notre Dame de Lourdes, MB R0G 1M0

T (204) 823-4037 | F (204) 248-2087

www.southernhealth.ca

TO/DESTINATAIRE : Acute/Home Care Directors & Clinical Managers, Pharmacy, Staff Development

CC : Lorraine Cassan, Regional Lead - Acute Care & Chief Nursing Officer
Dana Human, Regional Lead - Community & Continuing Care
Isaac Klassen, Manager - Physical Plant Services
Gerald Philippot, Clinical Manager -Emergency Medical Services

FROM/EXPÉDITEUR : Kerri Bolack, Manager Clinical Initiatives & Mentorship
Lori Wychopen, Site Manager-Carman Memorial Hospital
Tim Siran, Regional Manager – Supply Chain Services

DATE : January 30, 2023

PAGE(s) : 3

SUBJECT/OBJET : B Braun Infusomat Pump Implementation and Consumables

Southern Health-Santé Sud will begin implementation of the B Braun large volume infusion pumps in a staged approach across 13 facilities beginning the week of February 6th and continuing throughout the month. The training schedule, course registration instructions, Go-Live dates, and associated resources are available on the HPS [B Braun Implementation](#) page.

All primary infusion sets are proprietary to B Braun. Please review the chart below showing 'like for like' comparisons for consumables. Clinicians should note the clinical applications of each set and only order those that are applicable to your care area. The week your site is to go live, Materials Management will automatically send 4 weeks worth of the new infusion sets to you based on your previous usage. Beyond this initial supply, you should submit orders for your area as normal; noting the internal order number (SKUs) are not changing. All Baxter primary infusion sets will need to be removed from supply rooms and crash carts and exchanged for B Braun sets. Once you have gone live with the new tubing, please send the obsolete Baxter tubing back to Stores.











Baxter products that do not run through the pump itself should continue to be stocked; including secondary medication sets, extension sets, extension sets with filter, and needleless connectors, etc.







The Community Cancer Programs (BTHC, PDGH, and BRHC) will continue to use their current pumps/infusion sets until a Provincial Oncology Drug Library has been developed.

The SH-SS Parental IV Drug Monographs have also been updated to align with the launch of the new B Braun infusion pump and its onboard drug library. The Pharmacy Department has circulated print copies of the updated monographs to facilities for insertion into each local Parental Drug Monograph binder on the day of Go-Live. The [electronic version](#) of the monographs will also be updated on the HPS.

January 30, 2023

B Braun Infusion Set Conversion Chart

Current Product Mat Man #	Current Product Image	Item Description & B Braun Code	Future Product Image	Future Mat Man #	Clinical Use/Comments
00899 (Baxter 1C8109)		INFUSOMAT SPACE PUMP IV SET CA24 (363904)		00899	Straight set with no filter. Use during administration of all medications which you cannot "piggyback" into i.e. Heparin, Insulin, Amiodarone. If the medication requires a filter place an Extension set with 0.22 micron filter (Mat Man # 08955) at distal end.
00377 (Baxter 2C8419)		UNIV 15 DROP PUMP SET W/CARESITE CA24 ASV NON-DEHP (363421)		00377	2 ports with no filter. Can be used interchangeably for 00767
00767 (Baxter 2C8519)		UNIV 15 DROP PUMP SET 2 CARESITE CA24 ASV NON-DEHP (BB8250718SP)		00767	2 ports with no filter. Can be used interchangeably for 00377
00084 (Baxter 2C8537)		UNIV 15 DROP PUMP SET 3 CARESITE CA24 ASV NON-DEHP (363430)		00084	3 ports with no filter. Could be used to administer medications in medial port if mini bag not available or client fluid compromised
09283 (Baxter 2H8480)		UNIV 15 DR UV RESIST PUMP SETCA24 ASV NON-DEHP (363410)		09283	1 port with 0.2 micron filter. Used for Aminosyns & any medications that require a 0.2-micron filter as per drug monograph. If it is recommended

					that the medication should not be run with anything else (ie Amiodarone) use the Straight Set Mat Man # 00899 and an Extension 2H5660 non-DEHP 0.22 Micron Set, Mat Man # 08955
00378 (Baxter 2H8486)		UNIV 15 DR PUMP SET 1.2FIL 2CRSTE CA24 ASV NON-DEHP (363424)		00378	Used for IV Lipid Infusions. Best practice states these require a filter of at least 1.2 micron to prevent clogging.
00769 (Baxter 2H8864)		60 DROP METRISET PUMP 3 CARESITES CA24 ASV NON DEHP (363422)		00769	Buretrol set; may be used in pediatric care
01358 (Baxter JC7751)		Y-TYPE BLOOD PUMP SET CARESITE CA24 ASV NON-DEHP (363419)		01358	Administration of Blood

*ASV=Anti-Syphon Valve

ORDERING INFORMATION: Please remember you will all receive 4 weeks supply for Go-Live. Follow your usual ordering process for your site and update your ordering lists going forward.

If you have any questions related to clinical matters please contact Lori Wychopen at lwychopen@southernhealth.ca , Kerri Bolack at kbolack@southernhealth.ca or for purchasing/material management inquiries reach out to Tim Siran at tsiran@southernhealth.ca.