

Supply Chain Disruption Alert - BACKORDER

23 September 2024

SH-SS Revision 26 September 2024

This information is to provide supply updates for your awareness.

In an ongoing effort to mitigate supply disruption, Supply Chain Management-Shared Services (SCMSS) Value Analysis along with Clinical Change and Patient Safety Leads have identified the following replacement as **ALTERNATE**:



(15 drops/mL) Universal spike, roller clamp, free flow protection clamp, slide clamp, SPIN-LOCK® connector. Removable anti-siphon valve. Not made with DEHP or natural rubber latex. Priming volume: 22 mL, Length: 123 in. (approx. 312 cm)

ALTERNATE #1



(15 drops/mL) Universal spike, slide clamp, pressure limited check valve, two CARESITE luer access devices approx. 88 in. and 6 in. above distal end, free flow protection clamp, roller clamp, slide clamp, backcheck valve at distal Y-site, SPIN-LOCK connector. Removable anti-siphon valve. Not made with DEHP or natural rubber latex. Priming volume: 22 mL, Length: 120 in. (approx. 305 cm)

ALTERNATE #2



(15 drops/mL) Universal spike, slide clamp, pressure limited check valve, three CARESITE luer access devices approx. 104 in., 70 in. and 6 in. above distal end, free flow protection clamp, backcheck valve, roller clamp, SPIN-LOCK connector. Removable anti-siphon valve. Not made with DEHP or natural rubber latex. Priming volume: 23 mL, Length: 123 in. (approx. 312 cm)

**to be used in conjunction with: **

SecureSeal® Sterile Tamper-Evident Blue Foil IV Seals



BACKORDER

SAP 358165

SET IV STR 15DR/ML 22ML 123" ASV

Cardinal Product Code BB363904 Supplied 24EA/CA SH-SS Order # 00899

Replacement- Alternates **

SAP 328697

SET IV 2 INJ SITE 15DR/ML 22ML 120" ASV

Cardinal Product Code BB363421
Supplied 24EA/CA

<u>SH-SS Order</u> # 00377

-OR-

SAP 328694

SET IV 3 INJ SITE 15DR/ML 23ML 123" ASV

Cardinal product code BB363430
Supplied 24EA/CA

SH-SS Order # 00084

along with the following product to be clinically acceptable

SAP 285335

seal IV bag port blue

Healthmark product code WI-10-200 1000EA/BX

Blue Foil Port Seals is our recommended solution to cover the Care Ports in order to be Clinically Acceptable